

# ROACH User Manual

## *Tutorial*

### **Make a script for starting python**

Log into the roach machine with:

```
ssh -Y myusername@164.54.85.124
```

To start the roach software in a Linux box at APS you must run ipython, the Enthought Distribution supplied by APS. In your home linux directory make a script like this using the vi editor.

```
vi runipython
```

To enter lines into the file, hit *i*

Copy this line into the editor terminal:

```
/APShare/epd/rh6-x86_64/bin/ipython -pylab
```

Hit this key sequence to save, quit.

```
ESC
```

```
:w
```

```
:q
```

Make the script executable with

```
chmod +x runipython
```

### **Starting ROACH software in python**

Start python


```
./runipython
```

In the python screen, cd to your install of the roach software. For now we assume it is in

ROACH/projects. In the python terminal type the following lines:

```
cd ROACH
cd projects
execfile('natAnalGui.py')
main()
```


You should see this window open:



To get started, hit *Connect Roach* button. You may see some errors on the terminal, ignore them. They will get fixed eventually. Once you see *connction established* on the bottom of the window, you can search for resonators.

## Finding and making a list of resonators

Hit the *Sweep* Tab. You will see this:



Type in the name of your device, above it is M206. This will keep your data associated with your device name.


We will now find resonators on the device.

Set Span to 200MHz

Set Center Freq to 2500MHz

Hit *Start Sweep*

You will see plots of the resonators if they exist:



The plots above from top to bottom are

Amplitude in rms normalized to 1.0, Phase in radians, I, Q normalized to 1.0, rms. The x axis is frequency in Hz.


If you see no resonators,

Hit *Stop Sweep*

Increment Center(MHz) by hitting up arrow. Set to 2700MHz

Hit *Start Sweep*

Notice we see resonators below:




Let is ZOOM in.

Hit *Stop Sweep*

Set *Span* to 100MHz

Set *Center* to 2680MHz

Hit *Start Sweep*




Now we will make a list of resonators.

Hit *Settings* Tab.

Under *Plot Type*, select *MagPhase*

You should see:




Hit *Sweep* Tab


Hit *StopSweep*

Hit *extractRes*

A window comes up:



This shows that we found resonators correctly. The main window will look like this:



Note that resonators have red *x* marking them on the plot. Also we see a list of resonators to the right.

Notice that the names are NULL. This is a little bug...Under "Sweep and you shall Find" retype the name M206. Then the list updates. Sorry about the bug.

If you have problems finding resonators, with 'x' being in wrong places, you can delete resonators from the list, clear the list, and add resonators manually.

Clear the List:

Hit *ClearREsList*

To Delete a resonator:

In the list, select the one you wish to delete


Hit *DelRes*

To add resonator manually

Check *AddRes*

Click in the TOP plot, the magnitude, on the resonator you wish to add

Selected resonator BELOW




Now we have a list of resonators.

Let is save the list before we lose them!!!

Hit *ResData* Tab

Type in a path and name of a file and hit *Save*

You may see a “cannot stat” error message in your terminal, Ignore it. It is trying to make backups of your files. Once you save several times, you will have several backups in the ROACH/projects directory where you have your ROACH software.




## Resonator Power Sweep

Now you can do a power sweep. Follow these steps:

1. Set U6CnstOutAtt to some number of dB. You can leave it at 20. This attenuator will be held constant during the power sweep.
2. Set AttOutStart (U7) to a small attenuation, like 0Db.
3. Set AttOutEnd (U7) to a large attenuation like 20dB
4. Set AttOutStep to 2dB.
5. Set AttinStart (U28) to  $\text{AttOutEnd}(U7) + 1$ . IN this case set it to 21dB.

What we did is set up a power sweep. Attenuator U6 will be held constant at 20dB. Attenuator U7 will setp from 0dB to 20dB in 2dB steps. U28 will steo from 21dB to 1 dB in 2dB steps. Here is a diagraph of how these attenuators are connected:




Now we will start the power sweep:

1. In the list on the right, select boxes of resonators you wish to sweep. OR- Hit *Check All* button.
2. Check Multiprocess

3. Check Power Sweep
4. Check RunFits if you wish. It is slow, so you may skip it. You can run them later.
5. Check GetNoise. You will then get some noise data on each resonator. You can skip this, as it takes time. It is not the best measurement of noise anyway. There is a better way to be shown later.
6. Check runIQVelocity. This is essential if you get noise data.

Here is our setup:



Now hit

*RunChecked*

to start the power sweep.

It takes some time. Do something else.

When done, hit *Save*

The system makes backups of your sweep data called *ROACH/projects/powersweep\_backup.h5*

Also, your file will be saved automatically from time to time during the sweep. If the software crashes you can retrieve data from your file later, or from backup files called:

*ROACH/projects/powersweep\_backup.h5*

*ROACH/projects/backup\_N.h5*

The smallest N is the newest backup.

Again, *Roach/projects* could be a different directory. It is the home directory where your roach software are installed.

Once you get your data, see the section below on **Hdf Data Analysis**.

## Streaming Noise and Pulse data

### Noise Data

For best results, the resonators should have been power swepted at at least one attenuation and IQVelocity should have been run, to find the resonator centers. See above.

To stream noise data into a file, this is the best way to measure noise.

1. Hit *ResData* tab
2. Select resonators you wish to stream. For ROACH, the cleanest measurement is one res at a time because we have not developed high speed data transfer from ROACH. You can read out several resonators, but blocks of data are lost.
3. Hit *Stream* tab.
4. Hit *MKIDs->ROACH*. This gets the info from the selected resonators and sets up the ROACH board accordingly.
5. Type in a path/filename.
6. Type number of seconds you wish you stream data.
7. Make sure *PulseDet* is NOT checked.
8. Hit *StartStream*

You will end up with an hdf file on your computer with streamed data. See **HDF Data Analysis**.

### Pulse Data


To collect pulse data:

1. Do steps 1-6 above under NOISE DATA.
2. Check *PulseDet*
3. Set *PulseDet Nstd* to 5.
4. Hit *MeasMeans*. This takes a short noise measurement on resonators and gets mean values of noise and phase, magnitude values from each resonator.
5. Hit *ProgPulseDet*. This will program the pulse detector, and take a test measurement for .5 seconds to get pulse rate. AT *PulseDet Nstd* =5, you will see pulses.
6. Reduce the pulse rate by setting *PulseDet Nstd* to a larger number like 7. This is number of noise standard deviations.
7. Hit *ProgPulseDet* a few times to get count rate. A zero count rate may mean you never see pulses. A count rate of 80 pulses /sec is good. You will see something.
8. After setting up the pulse detector hit *StartStream*.

9. You should have an hdf file with pulses in it.

## Making Plots in the GUI

Once you have done a power sweep you will see resonators with many traces:




Select a trace as shown in the figure above.

Hit *ResonatorPlots*. You will see the following plots come up:

1. Mag, Phase, I, Q vs. Frequency. Sweep with noise included.
2. I vs Q raw (blue )and transrotated (red)
3. IQ velocity vs. Frequency. (sweep data)
4. Noise magnitude and phase plots.
5. Power spectrum of noise.

Plots are shown below.




Also when you hit *ResonatorPlots*, you will get info on the terminal on the traces

```

[screen 0: tcsh] tmadden@roach: ~
File Edit View Search Terminal Help
Attens U6 15.000000 U7 8.000000 U28 13.000000
FFT gain 2047
FFT base freq 0.000000
Num Noise Traces 20
Lut Sine amp 0.600000
Fits Ran 0 Fits Err 0
Rough Center Freq 2652.150284MHz
maxIQvel=0.000644,maxIQvel freq=2652131984.000000
Trace 97865559 Freq range 2650651880.000000, 1464.000000, 2653648688.000000
Carrier F0 2710000000.000000
Datalength 2048
Aprox Freqs:
[ 0. ]
Span 2.996808 MHz
Fitted Fr 0.000000 GHz
Fitted Q 1.000000
Fitted width 0.000000 MHz
Fitted Start 0.000000 GHz, End 0.000000 GHz
Attens U6 15.000000 U7 9.000000 U28 12.000000
FFT gain 2047
FFT base freq 0.000000
Num Noise Traces 20
Lut Sine amp 0.600000
Fits Ran 0 Fits Err 0
Rough Center Freq 2652.150284MHz
maxIQvel=0.000097,maxIQvel freq=2652158336.000000

```


If you select the MKID, and hit *ResonatorPlots* you will get plots of all sweeps.




All sweeps of MKID 0

## Sweeping and Live Plots

We can set up the system to sweep over and over a resonator like a network analyzer by setting the Span and Center Frequency on the *Sweep* Tab. See below.



Note that the resonator has too much power being sourced. Simply goto the Settings Screen and set the attenuators.








Additionally it is possible to loop back the system by clicking

- RF Loop- Loops back the RF signal at GHz size of Mixer.

- BB Loop- loops back Baseband, leaving mixer out of circuit.
- FPGA Loop- digital loop back inside the FPGA. No analog circuitry involved at all.

There are several plots that can be displayed live:

- four Plots- mag, phase, I, Q vs frequency for sweeps.
- Mag, Phase sweep
- I and Q sweep
- Polar
- I vs Q



There are more plots, but they do not work with the Sweeping network analyzer firmware.

## **HDF Data Analysis**

In ipython:

*hdfshell()*

Or if you are in the roach gui you can do this following:

1. On *ResData* tab, hit the *hdfShell* button
2. OR- Exit the ROACH gui to get python, by hitting the X at top of window.
3. OR- On *Settings* tab, hit *Exit2Console*. You will then have access to python.
4. In python type *hdfshell()*

*open /dir/dir/myhdf.h5*

Opens the hdf file.

*close*

closes the hdf file

*q*

Exits hdfshell() and saves hdf file.

*help*

*man*

print help screen

## **Navigating HDF5 files**

This starts a shell for navigating hdf files.

Use linux commands to navigate.

Opening hdf file

*ls -rlt*

Lists files on disk.

*cd /somedir/somdir etc*

cds in linux to dir.

*open backup\_1.h5*

opens hdf file backup\_1.h5

List contents

*ls*

*ls -ds*

will print dataset values

*ls dirname*

will print contents in that group

*cd groupname*

will set pwd to group

*pwd*

shows current group

*rm groupdataset*

moves group or data to the a group called /trash

*mv oldname newname*

renames or moves data in the hdf5 file

*cat mydataset*

*more dataset*

These commands show what is in a dataset.

## **Searching HDF files**

Use linux style find command:

*find* Finds everything at current dir.

*find .* Same as above

*find Device\_M206* finds everthing under dir Device\_M206

*find Device\_M206/Resonator\_1 -name fbase*

Finds everything under Device\_M206/Resonator\_1 with fbase in the name.

*find Device\_M206/Resonator\_1 -name fbase -ds*

Finds everything under Device\_M206/Resonator\_1 with fbase in the name, prints dataset values

*find Device\_M206/Resonator\_1 -name time -attr -ds*

Finds everyting in Device\_M206/Resonator\_1 with time in the name. Also searches attributes.  
Prints valuiies

*find Device\_M206/Resonator\_1 -name timefl -ds -out createtimefl*

Finds everyting in Device\_M206/Resonator\_1 with time in the name. Prints valuiies. These values will be stored into a new data set at current working dir, called createtimefl. You can use an absolute path for -out such as -out /mydir/mydir2/mydata

## **Plotting Datasets**

*plot xdata ydata*

*plot ydata*

plots dataset, 4k points

*plot -cs b xdata ydata*

plots x versys y data, 4k points

*plot -cs rx ydata*

Use python colorspec rx and plot ydata

*plot -ts ydata*

plot ydata events according to timestamp. Useful for pulses

*plot -l 40000 ydata*

plot 40000 points of ydata

*clf*

clears plot

## ***Making IQ Plots***

Because data may not be in the correct form for an IQ plot, the *pliq* command exists. Example, if we have mag and phase data, we must first convert to rectangular coordinates before making an IQ plot. Also, if data is stored in a multi-dimensional array, there is no easy way to extract the slabs of data from the dataset. For example, often resonator sweeps are stored as `iddata[2][L]` for L point sweeps.

*pliq -iq iqdata*

*pliq -mag magnitudedata -ph phasedata*

*pliq -i idata -q qdata*

For i and q data in separate

*pliq -c Chan+00000*

For a Chan\_XXXXX directory with magnitude and phase data.

*Pliq -l 40000 -i idata -q qdata*

Plot 40k points.

*Pliq -cs rx -iq iqdata*

plot iqdata is red x's

## **Noisie plots**

```
welch /Chan_XXXX/phase
```

makes PSD plot of phase noise. For looking at 1/f noise

## **Trans-Rotation of Resonator Sweep and Noise/Pulse data.**

To rotate data based on a sweep:

```
transrot -r /sweeps/ResData_4 -c /Chan_00000
```

Take sweep data in ResData\_4, use it to rotate the magnitude and phase in Chan\_00000. You must use the correct resonator sweep that corresponds to the Chan. To find out, match the rf frequency or baseband frequency of noise and sweep data.

Datasets magnitude\_tr and phase\_tr are created in the Chan\_XXXXX directory of transrotated data.

```
pltrot -r /sweeps/ResData_4 -c /Chan_00000
```

Plots resonator as a circle on IQ plane, plots noise cloud and pulses on the circle from Chan directory.

## **Making Notes**

To store a note in hdf file:

```
note mynote.txt
```

You then type a note. Hit ctrl-D to end note. The note is now stored.

To see note:

```
cat mynote.txt
```

## **Extracting Data to Python**

To extract a data set to python as a variable,

```
topy /mygroup/mydataset mypyvarname
```

The dataset is created as a python variable as a global, and in namespace myvars.

That is a global called mypyvarname is created and myvars. Mypyvarname is created.

```
py print mypyvarname
```

py will run a python command. Py and follow with any python command. Globals will be seen in the scope of the command. If you exit hdfshell() with 'q' then you are returned to python interpreter. The global py variables will still exist.

### Plotting data versus timestamps

Each resonator sweep has a timestamp when it was swept. The data sets are called

```
createtime
```

List of [year, month, day, hour, min, second]

```
createtimefl
```

float point number of seconds since epoch. See python time module for explanation. This is used as a timestamp when plotting data vs timestamp.

For some MKID, we can create a dataset of all sweep timestamps

```
find Device_M206/Resonator_9 -name createtimefl -ds -out Device_M206/Resonator_9/createtimefl
```

We find all create time stamps for that resonator for all sweep traces. These timestamps are put into dataset at

```
Device_M206/Resonator_9/createtimefl
```

You can plot this data:

```
plot Device_M206/Resonator_9/createtimefl
```

Or you can now get some other data point like the center freq of resonator, most useful is `maxIQvel_freq`. That is the res center freq based on max IQ velocity.

We make a new data set for centerfrequencies.

```
find Device_M206/Resonator_9 -name maxIQvel_freq -ds -out Device_M206/Resonator_9/centerfreqs
```

Now we have a centerfreqs dataset w/ center freq for each trace.

To plot centerfreq vs sweep time:

```
plot -cs x Device_M206/Resonator_9/createtimefl Device_M206/Resonator_9/centerfreq
```

## ***Examining Resonator Sweep Data***

Because taking sweeps and noise to characterize resonators creates much data, a special program called `resView` exists solely for dealing with resonator sweep data.

To start `resview`, resonator data must be loaded into a data structure internal to the python interpreter. In `hdfshell`, there are two types of resonator data:

1) Single sweep data from streamed noise

```
/sweeps/
```

2) Multi-Sweep data from sweeping resonators w/ power sweeps etc.

```
/Device_XXXX
```

*If there is a sweeps directory,*

```
cd sweeps
```

```
loadres
```

```
resview
```

Each resonator sweep will be loaded as a separate MKID, with one sweep each on `loadres`.

Resview starts the resonator viewer.

The structure of the sweeps director is as follows:

/sweeps/ResData\_0

/sweeps/ResData\_1

etc.

Each ResData\_N directory is a sweep of a separate MKID. IN this type of data file, each MKID has a single sweep.

*If there is a Device\_XXXX directory,*

*cd to root directory above Device\_XXXX dir.*

*loadres*

*resview*

The structure of Device\_XXXX is as follows:

Device\_XXXX/Resonator\_1

Device\_XXXX/Resonator\_2

etc

The Resonator\_X directories coorespond to indivitual MKIDs. Each MKID may hagne manuy sweeps.

Device\_XXXX/Resonator\_1/ResData\_0

Device\_XXXX/Resonator\_1/ResData\_1

etc.

The ResData\_N directories coorespond to each sweep trace for Resonaotr\_1.

After calling loadres,

*resview*

The resonator viewer program is started. You will see a prompt, ResView>>.

*q*

Exit resview back to hdfShell

## ***Resvview Commands***

Resview has the concept of current MKID and current trace. The idea is that we can step through MKIDs and traces and plot each one easily.

*q*

quit resview return to hdfshell

save

save hdf file.

help

man

List resview commands

list

lists MKIDs and all traces

<return>

Plot current trace

<space><return>

Increment to next trace and plot

+ <return>

Increment to next trace, plot, but do not clear plots

hold

nohold

Turn off automatic clear plots with hold. Nohold turns on auto clear of plots.

nmkid

goto Next MKID

st

Goto 1<sup>st</sup> MKID and 1<sup>st</sup> trace

plmkid

Plot all traces of this MKID. Need to call hold first so plot is not cleared after drawing.

fit

Run fits on current trace

clf

Clear plot, replot current trace

-<return>

Goto previous trace, plot

pr

Plot all fitting curves for this trace

iqvel

For this MKID, use all traces to make IQVelocity plot versus attenuation. This is for finding preferred attenuation for this MKID

prefat

Set preferred attenuation for this MKID. This is then stored in HDF File, and can be used later for taking data.

flatph

If phase is increasing too much w./ frequency due to wrong time delay in data, this sets time delay

properly for this trace. Makes fitting easier

fitprefat

fit resonator trace that is at preferred atten for this MKID. Finds correct trace based on atten, and fits.

Fitprefatall

Fit all MKIDs for traces at preferred atten.

Findat

find trace based on U7 atten

findat2

find trace for this MKID based on total atten, U6,7, and LUT amp.

plnoise

Turns on plotting of noise traces

nonoise

Turns off plotting of noise

## ***theory of Operation***

The ROACH runs two firmware designs:

1. Network Analyzer- A network analyzer design that can sweep resonators and read one resonator at a time.
2. FFT Analyzer- A code that uses FFTs for reading many channels at a time. Up to 256 channels can be read out.

The software interface for the ROACH is a collection of python scripts. Here is a list of the python scripts:

1. `t_brdconfig.py`- functions for IF board control, testing functions, network connection to ROACH, roach register setup. Control of anritsu freq. Generator.
2. `netAnaluzer.py`- control of network analyzer firmware. `NetAnaluzer` class is defined.
3. `fftanalyzer.py`, `fftanalyzerd.py`, `fftanalyzери.py`. These are codes for control of fft analyzer FW. They all inherit from `network analyzer` class. These three class are a line if inheritance. `Fftanaluzeri` inherits `fftanaluzerd`, which inherits `fftanaluzer`, which inherits `netanaluzer`. There were several versions of the fft formware, and each class represents major versions. The `i` version is the current version.
- 4 `fitters.py`- Resonator trace data objects, `MKID_list`, the current liost of MKIDs. `MKID opbjhects`, HDF read/write functions. `Fitter` class, which has all fitting functions. Several plotting functions.
- 5 `resciew,.py`- HDF file tool. A linux-shell that runs in ipython for navigatin HDF files.
- 6 `natAnaluzer.py`- GUI for controlling the ROPACH.
- 7, `roachMatlab.py`- a debygging script for allowing all py scripts to talk to matlab simulator instead of ROACH board. This is for testing opython codes with FPGA sumulation.
8. `controlScripts.py`, Misc. Py scripts for controlling roach. Many degugging scripts. Scripts for muylti-processing of fits. Sets up fitting queeus for multi process.
9. `roachpref.py`- preferences for roach gui. Edit this for your favorite gui setup.
- 10 `doit.py`- pyt anything you want in here. You can execute it from the gui.
- 11 `mpfit.py`- needed for `fitters.py`

Deprecated unused scripts:

`dumpView.py`

`fshift.py`

`sim.py`

`test.py`

`debug.py`

`ljh.py`

`doit.py`

`ftAnalyzere.py`

`runFits.py`

`test2.py`

## ***theory of operation of netanalyzer FW***

Registers

Data format returned

## ***theory of operation of FFT FW***

Registers

Data format in Fft output memory

## ***Streamed Data Readout***

## ***ROACH settings***

RF Switch settings

OSC settings

anritsu control

report()

Roach board ethernet settings

Sending FW to the roach board. /boffiles

Running Fw on roach.

Ioreg interface on roach

logging into roach

## ***Finding resonators***

**the Resonator List**

## ***Sweeping resonators***

## ***FFT Readout***

## ***Streaming data***

## ***Hdf files***

Dump files

## ***Debugging data- pckle***

## ***Useful Roach py scripts***

## **Controlling network analyzer from python**

## **Controlling FFTs from python**

Useful plots

## **Streams from python**

Resonator objects and fits

useful plots